PHILIPS

RESPIRONICS

Trilogy Evo

Quick start guide for caregivers

Overview		4
Your prescriptic	ons	5
Available circui	t options	6
Getting started	and the second se	10
Therapy views	And the second se	12
Prescription set	ttings	13
Options		15
How to custom	ize your background	16
Alarms		17
USB data		20
Device mainter	nance	22

Welcome to your Trilogy Evo ventilator from Philips Respironics. Trilogy Evo is designed for home, on-the-go, and hospital applications, for children and adults, so it can transition with patients through changes in their conditions or their environment.

Trilogy Evo features include:

- · Easy to use with a large (8") touch screen
- Up to 15 hours of power¹ for peace-of-mind during extended trips when away from home
- Less than 13 lbs (5.8 kg) including internal and integrated detachable battery, for the perfect blend of weight and power
- Increased durability and improved mobility with in-use carrying case, or wheel chair and roll-stand mounts²

So your ventilator can stay with you wherever you go.

- Run time mode settings for Trilogy Evo: A/C-VC mode ActivePAP circuit, PEEP of 3 cmH2O and Vt 800ml
- 2. Available for purchase

Adjusting to using your new ventilator, Trilogy Evo, will be easy.

- Settings to synchronize with your breathing pattern (inhale and exhale), for a natural and comfortable experience
- Multiple settings for timing, volume, and pressure, so it can be customized for your needs now and in the future
- Available support and training with additional information

This guide provides important information that will help you use the device. *Refer to the caregiver manual for complete details.*



Overview



			Circuit to us	a		Acc	cessories	
Prescription name	Passive see p.6	ActivePAP see p.7	Active Flow see p.8	Dual Limb see p.9	MPV see p.10	Humidifier	0 0	Other

Additional notes:

Overview

prescriptions:

Your Trilogy Evo has been programmed with _

Available circuit options Passive circuit



A. Connect the bacteria filter on the circuit to the inspiratory port.

Active PAP circuit



- A. Connect the bacteria filter on the circuit to the inspiratory port.
- B. Connect the proximal pressure line (wider diameter than active exhalation valve line) to the proximal pressure port.
- C. Connect the active exhalation valve pressure line to the active exhalation valve line connection.



Active Flow circuit



- A. Connect the bacteria filter on the circuit to the inspiratory port.
- **B.** Connect the proximal pressure line (wider diameter than active exhalation valve line) to the proximal pressure port.
- **C.** Connect the active exhalation valve pressure line to the active exhalation valve line connection.

- **D**. Attach the flow sensor cable to the flow sensor cable connector
- E. Attach the flow sensor to the active exhalation valve on the circuit



Dual Limb circuit



- A. Connect the bacteria filter end of the colored inspiration tube to the inspiratory port.
- **B.** Connect the proximal pressure line to the proximal pressure port.
- C. Install the active exhalation valve into the recessed AEV port. Press until both sides click into place.
- **D.** Attach the bacteria filter end of the clear expiration tube to the AEV.
- **E.** Attach the flow sensor cable to the flow sensor cable connector.
- **F.** Attach the flow sensor to the Y-shaped connector on the circuit.



Available circuit options (continued) MPV circuit



- A. Fully extend and straighten the circuit support arm. See the diagram below.
- B. Feed the circuit tube (15mm) through the center of the circuit support arm until it exits the other end.
- **C.** Attach the clamp to a wheelchair if required.
- D. Attach the reducer cuff and then the bacteria filter onto the device-end of the circuit tube.
- E. Connect the bacteria filter on the circuit to the inspiratory port on the Trilogy Evo.
- F. Attach the coupler and miniature flextube (optional) onto the circuit support arm before connecting your chosen patient interface.





Getting started Prepare for use



1. Plug in device

Use the AC cord provided to connect AC power. The green LED light next to the on/off (standby) button should be lit.



If AC power is not available, you can use battery or DC power.



2. Install filter

To install the air-inlet foam filter, pinch the filter as you press it into the filter cover as shown. Position it securely behind the top and bottom restraints.



3. Starting Trilogy Evo

Make sure all accessories, cords, and circuit connections are attached. Press the on/off (standby) button.

If starting on battery power, a notification will appear at the top of the screen along with a beep every 30 seconds. Tap the reset icon beside the notification to acknowledge.

Home standby window*

🚹 😳 d ¹ 2	Stanc Not Ventilat
Prescriptions	
A/C-PC AVAPS Passive Daytime	
A/C-PC AVAPS Passive Nighttime	
	tart Ventilation
	■ 2 ∧ 12:45

Your therapy prescriptions will be listed here, for selection

4. Select your prescription

5. Tap Start Ventilation to begin therapy

If you want to stop therapy, press the on/off (standby) button on the front panel, then tap Standby or Power Off.

Prescription menu

Views

Tap Views drop down arrow to access multiple monitoring options

Your prescription and circuit in use will be shown here. If your device is not ventilating, this area will read "Standby Not Ventilating"



When the current breath is triggered by the patient, this indicator appears filled (dark green)

Manometer

This bar will move based on the pressure being delivered

Battery time

Estimated remaining battery time, based on prescription in use and breathing pattern



Touchscreen lock

To prevent accidental therapy changes, use Touchscreen Lock, Lock the screen any time with the Status Bar shortcut shown here

Prescription settings Prescriptions settings window

Tap the settings icon to view more details about your prescriptions

Daytime	\sim	Start Ve	entilation		• Alarm
Circuit Passive Mode A/C-PC AVAPS Advanced	Tidal Volume 400 mL Breath Rate 15 BPM AVAPS Speed 5 cmH20/min	PC Min/Max 10/20 cmH20 Trigger Type Auto Trak	PEEP 5 cmH2O Trigger Sens. Auto	Insp. Time 1.5 s Rise Time 2	Switch tabs to view alar settings
Circuit			1	Using Default Calibration	
Type Passive	Size Adult	: (19-22mm)	Active Hum Off	idification	

Tap Circuit, Mode, or Advanced options on the left to see more information

Alarm settings

li 🌻 di	3				Not V	Standby entilating
Daytime	\sim	Start Ve	ntilation	© 🗘		
Circuit Passive	Tidal Volume Off/Off mL	MinVent 3.5/Off L/min	$\frac{\text{Resp. Rate}}{Off/45}\text{BPM}$	Circuit Disconnect 10 s		
A/C-PC AVAPS						
Advanced						
				e	^	12:45pm

Prescription settings (continued) Selecting a different prescription while ventilating

Option 1.

In the home window, tap the prescription drop-down arrow to access the prescription menu. Then, select a prescription to switch therapy



Option 2.

From the settings window, tap the drop-down arrow, select your new prescription, and then tap Switch Therapy

🔂 🔅 ය්ප			A/C-PC A	/APS Passive Daytime 〜	,
Daytime	Start Ve	ntilation	۵ ۵		
A/C-PC AVAPS Passive Daytime	me PC Min/Max - 10/20 cmH20	PEEP 5 cmH20	Insp. Time 1.5 s		
A/C-PC AVAPS Passive Nighttime	te Trigger Type Auto Trak	Trigger Sens. Auto	Rise Time 2		
Advanced AVAPS	Speed 20/min				
Circuit		Ustr	ng Default Calibration		
Type Passive	^{Size} Adult (19-22mm)	Active Humidi	fication		
			- 2	∧ 12:45pr	n

Circuit note: The circuit settings must be the same as the current prescription. If the circuit settings differ, place the device into standby to change the physical circuit. Then, select the prescription from the home screen to start ventilation.

Options Options window

Tap the options icon for the Options menu window

☆ 🖗				Not V	Standby
Options					
Device Options	>	Calibration & Setup	>		
Data Transfer	>	Alarm & Event Log	>		
Information	>	Prescription Preferences	>		
		Ī			
			• Z	^	12:45pm
		Within this window, options, run calibrati and view and work v	change d ions and vith data	evice tests,	

Device options

-	3	Off
Automatic Touchscreen Lock Screen Saver Off Breath	Date and Format MM/DD/YYYY	Time and Format 12-Hour
Bluetooth Not Enabled		

You are able to tap any of these device options and change if needed, including light bar, touchscreen lock, screen saver, and more



Alarms Disconnected Circuit example

This is an example of an alarm situation that may occur and how to respond.



Alarm silence button

Visually monitor the patient and ventilator at all times during an alarm silence period. Allowing alarms to continue without intervention may result in harm to the patient. If you want to silence the audible alarm for two minutes, press the Alarm Silence button

Resetting an alarm

Reconnect the circuit. When the alarm condition resolves, the alarm message will turn gray



Alarms (continued)

lcon	Description	Light and Sound Indicators
	High priority alarm	Light bar flashes red
		Audible alarm repeats rapidly
	Medium priority alarm	Light bar flashes yellow
	Low priority alarm	Audible alarm repeats moderately
	System message	Single beep
	Resolved alarm or system message	None

Alarm name	What happened?	What to do
Ventilator Inoperative	The system self-test indicates a failure or malfunction of a component. The failure causes therapy to stop or not meet essential performance criteria.	Provide an alternate method of ventilation then contact customer service.
Ventilator Service Required	This alarm occurs when the device cannot perform to specification, a backup safety feature is compromised, or the delivery of therapy is compromised. The device continues to function.	Contact customer service.
Obstruction	The ventilator detects an obstruction in the patient's inhalation path, exhalation path, or external flow sensor. The ventilator detects that the leak device is missing.	Check the circuit. Is it kinked or pinched? Check the bacterial filter. Is it blocked? If using an Active Flow or Dual Limb circuit, check the HME. Is it blocked? Is the leak device blocked or missing? Is the external flow sensor blocked?
Apnea	Time between patient- initiated breaths is more than the alarm setting.	 Is the circuit connected to the patient? Is there a leak or disconnect? Is the circuit kinked or pinched?
High Expiratory Pressure	During the expiratory phase, the delivered pressure exceeds the target patient pressure by 5cm H ₂ O or more.	Check the circuit. Is it kinked or pinched? Is the leak device blocked or occluded?

Alarm name	What happened?	What to do
High Inspiratory Pressure	During the inspiratory phase, the delivered pressure exceeds the target patient pressure by 5cm H ₂ O or more.	 Check the patient: Is the patient coughing or having excessive secretions? Is the patient having bronchospasms? Is the tracheotomy tube stable? Check the ventilator: Is the circuit kinked or pinched? Is the leak device or the exhalation device blocked? Are secretions in the HME?
Active Exhalation Valve Failed	The active exhalation valve is stuck closed.	 Single-limb active circuit: check all connections to the valve and check that the valve is clear. Dual limb circuit: check that the valve is clear. Replace the circuit or valve if necessary.
Check AEV Pilot Line	The active-exhalation pressure control line is not connected, becomes disconnected, or contains water droplets that affect the active exhalation valve line pressure reading.	Inspect the line. Empty or replace it if necessary.
Proximal Pressure Line Disconnected	The proximal pressure line is not connected.	 Is the line connected at both ends? Is the line clean and not tangled? Be sure the main circuit is connected and does not have large leaks Be sure the exhalation valve is intact.
Circuit Disconnected (MPV)	When using mouthpiece ventilation, the system detects the circuit is disconnected.	Reconnect the circuit.
Circuit Disconnected	The patient is not connected to the ventilator breathing circuit or there is a large leak.	 Is the circuit connected to the patient? Is the circuit connected to the ventilator? Does a large unplanned leak exist?
Low Minute Ventilation	The patient's minute ventilation is less than or equal to the alarm setting. Or, no breath has occurred for 15 seconds.	Is the circuit kinked or pinched? Does the circuit have a leak or disconnect? Remove excessive water from the tubing. Is the bacterial filter blocked or not connected? Is the leak device blocked or not connected? Check the patient

Alarms (continued)

Alarm name	What happened?	What to do
Low Respiratory Rate	The patient's respiratory rate is less than or equal to the low respiratory rate alarm setting. Or, no breath has occurred for 15 seconds.	 Is the circuit kinked or pinched? Does the circuit have a leak? Is the circuit connected? Check the patient.
High Inspiratory Pressure	Applies to volume modes. The measured patient pressure exceeds the High Inspiratory Pressure setting.	 Check the patient: Is the patient coughing or has excessive secretions? Is the patient having bronchospasms? Is the tracheotomy tube stable? Check the ventilator: Is the circuit kinked or pinched? Is the leak device or the exhalation device blocked? Are secretions in the HME?
No Trigger	Applies to mouthpiece ventilation. The time between patient-triggered breaths is greater than the No Trigger alarm setting.	 Check the patient Is the circuit within reach? Is the circuit intact? Is the PEEP high enough to generate sufficient signal flow when the mouthpiece is touched?

USB data

If your homecare provider or clinician asks you to download your ventilator data to USB in order to review your therapy, carefully follow the next steps:





- Connect the USB drive into any of the 2 USB connectors on either sides of the device.
- 2. In the menu bar, tap the CP options icon and then tap Data Transfer.
- Tap Export Data–USB to start the transfer to the USB drive.
- 4. Tap USB in the pop-up to confirm the destination for the export.
- 5. Tap Yes in the pop-up to confirm export ventilator data to USB.
- 6. Wait for Export Data-USB to display Success before removing USB stick.

Notes

Device maintenance* Cleaning the device

Frequency

Clean Trilogy Evo's exterior surface weekly

Supplies

- 1 lint-free cloth
- 1 soft-bristle brush
- Liquid dishwashing detergent solution: 1 teaspoon of liquid dishwashing detergent per gallon of warm water



Instructions



 To clean the exterior, first turn the device off and disconnect it from the power source detaching all accessories and connectors.



 Use a dampened lint-free cloth (not dripping) with detergent solution to clean the exterior of the enclosure.



 Use a soft-bristle brush in the areas around the screen, buttons, and any other areas where soil may be difficult to remove. Ensure you remove all visible soil.



 Use a dampened lint-free cloth (not dripping) with clear water to remove all detergent residue.



- 5. Use a dry lint-free cloth to dry the enclosure.
- Inspect the device for cleanliness and repeat the cleaning steps until the surfaces are visibly clean.



 Inspect the device for damage after cleaning. If any parts are damaged, contact customer service.

*For accessories cleaning/disinfection instructions vary based on the type of accessories you use. Please refer to your IFU for how to clean/dispose of your circuits.

Cleaning the battery



- 1. Open the detachable battery access door.
- Lift the battery handle and pull it to remove it from the bay and follow the instructions for cleaning the device (pg. 22).
- 3. Allow the battery to air dry completely and then inspect it for damage after cleaning.
- Slide the battery back into the bay until you hear a click.

Rinsing the air-inlet foam filter

The air-inlet foam filter is the gray foam located on the back panel. It protects Trilogy Evo from dirt and dust. This filter is for single patient use. Only use Philips Respironicssupplied filters. Ventilation can continue while you are replacing the filter.

Frequency

Rinse monthly and replace every six months.





Instructions

- 1. Ensure you have a replacement filter nearby.
- 2. Pinch the filter and pull it out of the filter cover.
- Insert the clean replacement filter into the filter cover. Ensure it is positioned securely.
- Visually inspect the filter you just removed from the device – if it is damaged, discard it. Otherwise, rinse the dirty filter in clear water. Inspect the filter for cleanliness and repeat until the filter is clean.
- 5. Allow the filter to air dry completely before reinstalling it.

Notes

Reorder information

Description





Respironics Inc. 1001 Murry Ridge Lane Murrysville, PA 15668 USA

Respironics Deutschland GmbH & Co. KG Gewerbestrasse 17 82211 Herrsching, Germany



REF 1141755 1141755 ROO ZL 05/06/2019

EC REP